

# ESTELITE <sup>Sigma</sup> QUICK

Universal Restorative Resin Composite

Estelite® Sigma Quick is a universal, 100% spherically filled composite, with a filler weight of 82% (71% volume). With Estelite Sigma Quick the best just got better; offering all the same characteristics of Estelite Sigma with the addition of **three new benefits**.

Compared to Estelite Sigma, Estelite Sigma Quick offers a reduced curing time (by one third) and an extended working time (90 seconds) under the dental light by utilizing Tokuyama's patented innovative initiator system "Radical Amplified Photopolymerization (RAP) Technology". Estelite Sigma Quick also has three new shades; A5, OPA2 (Opaque), and WE (White Enamel), a new bleach shade.




## Benefits


- Extended working time (90 sec.)
- Quick (reduced) curing time
- Outstanding polishability
- High gloss retention over time
- Wide shade matching range/"Chameleon effect"
- High wear resistance
- Less wear to opposing teeth
- Excellent radiopacity
- Minimal shrinkage
- Available in 20 shades

## Indications

- Anterior and posterior restorations
- Composite veneer
- Diastema closure
- Composite/porcelain repair



Estelite® Sigma Quick



Filtek™ Supreme\*



Tetric® EvoCeram\*\*




Venus™\*\*


\*Not a registered trademark of Tokuyama Dental Corp.

\*\*Courtesy of Tokuyama Dental Corporation, Tsukuba Research Laboratory, Japan


Working Time (sec)\*\*



Polymerization Shrinkage (vol%)\*\*



Glossiness (%)\*\*



# ESTELITE $\Sigma$ QUICK

Universal Restorative Resin Composite

## Products Available

### Estelite® Sigma Quick Syringe

● 3.8g/syringe



### Estelite® Sigma Quick PLT (Pre-Loaded Tip)

● 0.2g/tip x 20




### Estelite® Sigma Quick PLT Mini Kit

● 20 PLTs (A1, A2, A3 & OA2) 0.2g/tip x 5



### Estelite® Sigma Quick Shade Guide



### SPECIALTY SHADES



## Estelite® Sigma Quick Products

ITEM #	ITEM DESCRIPTION	ITEM #	ITEM DESCRIPTION
13186	Estelite® Sigma Quick Shade Guide	13381	Estelite® Sigma Quick PLT Mini Kit
13110	Estelite® Sigma Quick Syringe (A1)	13310	Estelite® Sigma Quick PLT (A1)
13111	Estelite® Sigma Quick Syringe (A2)	13311	Estelite® Sigma Quick PLT (A2)
13112	Estelite® Sigma Quick Syringe (A3)	13312	Estelite® Sigma Quick PLT (A3)
13113	Estelite® Sigma Quick Syringe (A3.5)	13313	Estelite® Sigma Quick PLT (A3.5)
13114	Estelite® Sigma Quick Syringe (A4)	13314	Estelite® Sigma Quick PLT (A4)
13115	Estelite® Sigma Quick Syringe (A5)	13315	Estelite® Sigma Quick PLT (A5)
13116	Estelite® Sigma Quick Syringe (B1)	13316	Estelite® Sigma Quick PLT (B1)
13117	Estelite® Sigma Quick Syringe (B2)	13317	Estelite® Sigma Quick PLT (B2)
13118	Estelite® Sigma Quick Syringe (B3)	13318	Estelite® Sigma Quick PLT (B3)
13119	Estelite® Sigma Quick Syringe (B4)	13319	Estelite® Sigma Quick PLT (B4)
13120	Estelite® Sigma Quick Syringe (C1)	13320	Estelite® Sigma Quick PLT (C1)
13121	Estelite® Sigma Quick Syringe (C2)	13321	Estelite® Sigma Quick PLT (C2)
13122	Estelite® Sigma Quick Syringe (C3)	13322	Estelite® Sigma Quick PLT (C3)
13123	Estelite® Sigma Quick Syringe (OA1)	13323	Estelite® Sigma Quick PLT (OA1)
13124	Estelite® Sigma Quick Syringe (OA2)	13324	Estelite® Sigma Quick PLT (OA2)
13125	Estelite® Sigma Quick Syringe (OA3)	13325	Estelite® Sigma Quick PLT (OA3)
13126	Estelite® Sigma Quick Syringe (OPA2)	13326	Estelite® Sigma Quick PLT (OPA2)
13127	Estelite® Sigma Quick Syringe (CE)	13327	Estelite® Sigma Quick PLT (CE)
13128	Estelite® Sigma Quick Syringe (BW)	13328	Estelite® Sigma Quick PLT (BW)
13129	Estelite® Sigma Quick Syringe (WE)	13329	Estelite® Sigma Quick PLT (WE)

For more information, please call 1.877.378.3548 or visit [www.tokuyama-us.com](http://www.tokuyama-us.com)

