



Estelite Omega... Redefining Esthetic Restorations

Estelite Omega is the newest in luxury and innovation, created for the artist within you. **Estelite Omega** is an elite universal, 100% spherically filled supra-nano composite with a filler weight of 82% (78% volume). Conceptualized and crafted to be specifically used with the next stage in esthetic restorations, **Estelite Omega** is the ultimate composite system to perform the latest in layering techniques. Developed in collaboration with Dr. Newton Fahl, Jr., **Estelite Omega** comes to you with praise amongst leading and esteemed doctors in the dental industry. With only 11 shades, your boundaries to perfect your masterpiece are limitless. Finally, a brilliant composite to showcase your creativity and elevate your art!

Benefits:


- Available in 11 Shades "Do More with Less"
- Optimal Esthetics
- Outstanding Polishability
- Exceptional Gloss Retention
- 90 Sec. Working Time
- Quick Curing Time
- High Flexural and Compressive Strength

Indications:


- Anterior and Posterior Restorations
- Composite Veneer
- Diastema Closure
- Composite/Porcelain Repair



Estelite Omega



Filtek Supreme Plus



*Courtesy of Tokuyama Dental Corporation, Tsukuba Research Laboratory, Japan

Estelite Omega Gallery

Clinical Cases as Performed by Dr. Newton Fahl, Jr.

Class IV Restoration




Direct Veneer Restoration



Class V Restoration



Diastema Closure



Palette of Shades

Estelite Omega is now available in only 11 shades... in hues and variances, providing you with the tools you need for natural and life-like restorations. Now you really can do more with less! Value simplicity at its best, without compromising esthetics.



Customize Your Palette

Estelite Omega's Custom Shade Guide allows you to build and personalize according to your own preferences. No more guessing when choosing your composite shade with this valuable tool. Customize your palette and create the exact shade of your composite.



Indulge in Colors

Compliment your restorations with **Estelite Color**, used to add character to your masterpiece, perfecting the overall esthetics. **Estelite Color** takes you to new heights in creativity, which only the inner artist would appreciate. Available in 13 colors, **Estelite Color** will advance your creative touch.



ESTELITE OMEGA

Estelite Omega Product Availability



SYRINGE KITS

10701 Estelite Omega Syringe Deluxe Kit
 11 Syringes (4g) in 11 different shades:
 Enamel: EA1, EA2, EA3, EB1
 Effects: TRANS, MW
 Dentin: DA1, DA2, DA3
 Bleach: BL1, BL2
 Custom Shade Guide
 Artist's Brush (Double-Ended)
 DVD/CD (Clinical Cases & Clinical Guide)
 Restorative Technique Card
 Instructions



10702 Estelite Omega Syringe Basic Kit
 6 Syringes (4g) in 6 different shades:
 Enamel: EA1, EA2, EB1
 Effects: TRANS, MW
 Dentin: DA2
 Custom Shade Guide
 Artist's Brush (Double-Ended)
 DVD/CD (Clinical Cases & Clinical Guide)
 Restorative Technique Card
 Instructions



PLT KITS (Pre-Loaded Tips)

10711 Estelite Omega PLT Deluxe Kit
 100 PLTs (0.25g) in 11 different shades:
 10 ea Enamel: EA1, EA2, EA3, EB1
 5 ea Effects: TRANS, MW
 10 ea Dentin: DA1, DA2, DA3
 10 ea Bleach: BL1, BL2
 Custom Shade Guide
 Artist's Brush (Double-Ended)
 DVD/CD (Clinical Cases & Clinical Guide)
 Restorative Technique Card
 Instructions



10712 Estelite Omega PLT Basic Kit
 50 PLTs (0.25g) in 6 different shades:
 10 ea Enamel: EA1, EA2, EB1
 5 ea Effects: TRANS, MW
 10 ea Dentin: DA2
 Custom Shade Guide
 Artist's Brush (Double-Ended)
 DVD/CD (Clinical Cases & Clinical Guide)
 Restorative Technique Card
 Instructions



SYRINGE REFILLS

10721 Syringe EA1 (4g)
 10722 Syringe EA2 (4g)
 10723 Syringe EA3 (4g)
 10724 Syringe EB1 (4g)
 10725 Syringe TRANS (4g)
 10726 Syringe MW (4g)
 10727 Syringe DA1 (4g)
 10728 Syringe DA2 (4g)
 10729 Syringe DA3 (4g)
 10730 Syringe BL1 (4g)
 10731 Syringe BL2 (4g)



PLT REFILLS (Pre-Loaded Tips)

10741 PLT EA1 (20 x 0.25g)
 10742 PLT EA2 (20 x 0.25g)
 10743 PLT EA3 (20 x 0.25g)
 10744 PLT EB1 (20 x 0.25g)
 10745 PLT TRANS (20 x 0.25g)
 10746 PLT MW (20 x 0.25g)
 10747 PLT DA1 (20 x 0.25g)
 10748 PLT DA2 (20 x 0.25g)
 10749 PLT DA3 (20 x 0.25g)
 10750 PLT BL1 (20 x 0.25g)
 10751 PLT BL2 (20 x 0.25g)

ACCESSORIES


10762 Brush no. 24 (5 Brushes)
 10763 DVD/CD (Clinical Cases & Clinical Guide)

CUSTOM SHADE GUIDE

10761 Custom Shade Guide

ESTELITE COLOR

10901 Estelite Color Kit (13 x 0.9g)
 10911 Clear (0.9g)
 10912 White (0.9g)
 10913 Blue (0.9g)
 10914 Yellow (0.9g)
 10915 Gray (0.9g)
 10916 Ochre (0.9g)
 10917 Dark Brown (0.9g)
 10918 Red (0.9g)
 10919 Lavender (0.9g)
 10920 Pink Opaque (0.9g)
 10921 High Chroma Opaque (0.9g)
 10922 Med Chroma Opaque (0.9g)
 10923 Low Chroma Opaque (0.9g)



Tokuyama Dental America Inc.

740 Garden View Ct., Ste. 200
 Encinitas, CA 92024
 Toll Free: 1-877-ESTELIT(E) or 1-877-378-3548
 Phone: (760) 942-7211
 Fax: (760)942-7212
 Email: info@tokuyama-us.com
 www.esteliteomega.com

